

GPR91 Rabbit pAb

ABZ1304

Product Overview

Name	GPR91 Rabbit pAb
Accession(Primary)	Q9BXA5
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Pig
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:200-1:1000

Target information (Q9BXA5)

Gene Name SUCNR1

About Gene Symbol



