

FSP1/AIFM2 Rabbit pAb


ABZ1305

Product Overview

Name	FSP1/AIFM2 Rabbit pAb
Accession(Primary)	Q9BRQ8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:40000

Target information (Q9BRQ8)

Gene Name	AIFM2
About Gene Symbol	



Western blot analysis of FSP1 in COS7 lysates
using FSP1 antibody.