

MyD88 Rabbit pAb


ABZ1306

Product Overview

Name	MyD88 Rabbit pAb
Accession(Primary)	Q99836
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Pig
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; ICC-1:-1:200-1:1000; IF-1:50-1:200; ELISA-1:20000

Target information (Q99836)

Gene Name	MYD88
About Gene Symbol	MYD88 (MYD88 Innate Immune Signal Transduction Adaptor) is an innate immune sensor involved in innate immune signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellu



Western blot analysis of MyD88 in CHO-K1, Jurkat, 3T3, C6 and Hela lysates using MyD88 antibody.