

CMKLR1 Rabbit pAb


ABZ1307

Product Overview


Name	CMKLR1 Rabbit pAb
Accession(Primary)	Q99788
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (Q99788)


Gene Name	CMKLR1
About Gene Symbol	




Western blot analysis of CMKLR1 in various lysates using CMKLR1 antibody.




Western blot analysis of CMKLR1 in COLO205 lysates using ChemR23 antibody.



Western blot analysis of CMKLR1 in 293T lysates using ChemR23 antibody.



Western blot analysis of CMKLR1 in 293T lysates using CMKLR1 antibody.



Western blot analysis of CMKLR1 in HeLa and HUVEC lysates using CMKLR1 antibody. The lane on the right is blocked with the synthesized peptide.