

Collagen XII alpha 1 Rabbit pAb

ABZ1308


Product Overview

Name	Collagen XII alpha 1 Rabbit pAb
Accession(Primary)	Q99715
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:100-1:300; IF-1:200-1:1000

Target information (Q99715)

Gene Name COL12A1

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human muscle using Collagen XII alpha 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.