


EYA1/4 Rabbit pAb

ABZ1309


Product Overview

Name	EYA1/4 Rabbit pAb
Accession(Primary)	Q99502
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200


Target information (Q99502)

Gene Name EYA1/4**About Gene Symbol**


Immunohistochemistry analysis of paraffin-embedded Human kidney using EYA1/4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using EYA1/4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using EYA1/4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using EYA1/4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.