


ZFP91 Rabbit pAb

ABZ1311

Product Overview

Name	ZFP91 Rabbit pAb
Accession(Primary)	Q96JP5
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q96JP5)

Gene Name ZFP91**About Gene Symbol**

Western blot analysis of ZFP91 in hepg2, NT28 lysates using ZFP91 antibody.

