

VSIG2 Rabbit pAb

ABZ1312


Product Overview

Name	VSIG2 Rabbit pAb
Accession(Primary)	Q96IQ7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q96IQ7)

Gene Name VSIG2

About Gene Symbol



Western blot analysis of VSIG2 in hela lysates using VSIG2 antibody.

