

SIRT1 Rabbit pAb


ABZ1313

Product Overview

Name	SIRT1 Rabbit pAb
Accession(Primary)	Q96EB6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 82 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (Q96EB6)

Gene Name	SIRT1
About Gene Symbol	



Western blot analysis of SIRT1 in various lysates using SIRT1 antibody.