

PLD4 Rabbit pAb

ABZ1314


Product Overview

Name	PLD4 Rabbit pAb
Accession(Primary)	Q96BZ4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q96BZ4)

Gene Name PLD4

About Gene Symbol



1 customer's

2 mouse-brain

3 mouse-kidney

Western blot analysis of PLD4 in various lysates using PLD4 antibody.