

FGF11 Rabbit pAb

ABZ1315


Product Overview

Name	FGF11 Rabbit pAb
Accession(Primary)	Q92914
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q92914)

Gene Name FGF11

About Gene Symbol



Western blot analysis of FGF11 in NIH3T3, HeLa lysates using FGF11 antibody.

