


CCL4/MIP-1 beta Rabbit pAb


ABZ1317

Product Overview

Name	CCL4/MIP-1 beta Rabbit pAb
Accession(Primary)	Q8NHW4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 10 kDa; Observed MW: 10 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q8NHW4)

Gene Name CCL4L1**About Gene Symbol**



Western blot analysis of CCL4L in K562 lysates using MIP1b antibody.