

BEST2 Rabbit pAb

ABZ1318

Product Overview

Name	BEST2 Rabbit pAb
Accession(Primary)	Q8NFU1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:50-1:200

Target information (Q8NFU1)

Gene Name BEST2

About Gene Symbol

