

IL-17 RC Rabbit pAb

ABZ1319


Product Overview


Name	IL-17 RC Rabbit pAb
Accession(Primary)	Q8NAC3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 86 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q8NAC3)


Gene Name IL17RC

About Gene Symbol





Western blot analysis of IL-17 Receptor C in HuvEc lysates using IL-17RC antibody.



Western blot analysis of IL-17 Receptor C in RAW264.7 lysates using IL-17RC antibody.