


cGAS Rabbit pAb

ABZ1320

Product Overview

Name	cGAS Rabbit pAb
Accession(Primary)	Q8N884
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:500-1:2000

Target information (Q8N884)

Gene Name CGAS**About Gene Symbol**

Western blot analysis of cGAS in HT-29 and Lovo
lysates using cGAS antibody

