

FBXO16 Rabbit pAb

ABZ1322


Product Overview


Name	FBXO16 Rabbit pAb
Accession(Primary)	Q8IX29
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:200

Target information (Q8IX29)

Gene Name FBXO16**About Gene Symbol**

kDa





Immunohistochemistry analysis of paraffin-embedded Human liver cancer using FBXO16 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.