


FUNDC1 Rabbit pAb

ABZ1323


Product Overview

Name	FUNDC1 Rabbit pAb
Accession(Primary)	Q8IVP5
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:100-1:300

Target information (Q8IVP5)

Gene Name FUNDC1**About Gene Symbol**

Western blot analysis of FUNDC1 in AD293 lysates using FUNDC1 antibody.



Western blot analysis of FUNDC1 in AD293 lysates using FUNDC1 antibody.