

STING Rabbit pAb


ABZ1324

Product Overview

Name	STING Rabbit pAb
Accession(Primary)	Q86WV6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q86WV6)

Gene Name	STING1
About Gene Symbol	STING1 (Stimulator Of Interferon Response CGAMP Interactor 1) is an innate immune sensor involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely



Western blot analysis of Transmembrane Protein
173 in 3T3 lysates using STING antibody.