

MACC1 Rabbit pAb

ABZ1326

Product Overview

Name	MACC1 Rabbit pAb
Accession(Primary)	Q6ZN28
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 97 kDa; Observed MW: 97 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q6ZN28)

Gene Name MACC1

About Gene Symbol



Western blot analysis of MACC1 in L929 lysates using MACC1 antibody.

