

MG53 Rabbit pAb

ABZ1327

Product Overview

Name	MG53 Rabbit pAb
Accession(Primary)	Q6ZMU5
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:50-1:200


Target information (Q6ZMU5)

Gene Name


TRIM72

About Gene Symbol

94KD
66KD
45KD
35KD
26KD



Western blot analysis of MG53 in mouse skeletal muscle lysates using MG53 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse Pancreas Tissue using MG53 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.