

CD300LD Rabbit pAb

ABZ1328


Product Overview

Name	CD300LD Rabbit pAb
Accession(Primary)	Q6UXZ3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000


Target information (Q6UXZ3)

Gene Name CD300LD

About Gene Symbol



Western blot analysis of CD300d in NIH3T3
lysates using CD300d antibody.



Western blot analysis of CD300d in 3T3 lysates
using CD300d antibody.