

RFWD3 Rabbit pAb


ABZ1330

Product Overview


Name	RFWD3 Rabbit pAb
Accession(Primary)	Q6PCD5
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (Q6PCD5)


Gene Name	RFWD3
About Gene Symbol	



Western blot analysis of RFWD3 in 823 lysates using RFWD3 antibody.



Western blot analysis of RFWD3 in MCF-7 lysates using RFWD3 antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of RFWD3 in HepG2 lysates using RFWD3 antibody.