

GDF6 Rabbit pAb

ABZ1331


Product Overview

Name	GDF6 Rabbit pAb
Accession(Primary)	Q6KF10
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000


Target information (Q6KF10)

Gene Name GDF6


About Gene Symbol



Western blot analysis of GDF6 in rat kidney lysates using GDF6 antibody.



Western blot analysis of GDF6 in rat kidney lysates using GDF6 antibody.



Western blot analysis of GDF6 in 293T lysates using GDF6 antibody.