

TRIM68 Rabbit pAb


ABZ1333

Product Overview


Name	TRIM68 Rabbit pAb
Accession(Primary)	Q6AZZ1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (Q6AZZ1)

Gene Name	TRIM68
About Gene Symbol	



Western blot analysis of TRIM68 in Jurkat, K562 lysates using SS56 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using TRIM68 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.