

GORAB Rabbit pAb

ABZ1335

Product Overview

Name	GORAB Rabbit pAb
Accession(Primary)	Q5T7V8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q5T7V8)

Gene Name GORAB

About Gene Symbol



Western blot analysis of GORAB in AD293 lysates using GORAB antibody.

