

ASAH3L Rabbit pAb


ABZ1336

Product Overview


Name	ASAH3L Rabbit pAb
Accession(Primary)	Q5QJU3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (Q5QJU3)


Gene Name	ACER2
About Gene Symbol	




Western blot analysis of ASAHL in COS7 lysates using ASAHL antibody.



Immunohistochemistry analysis of paraffin-embedded Human breastcancer using ASAHL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human breastcancer using ASAHL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human liver using ASAH3L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.