

TNFAIP8L3 Rabbit pAb

ABZ1337


Product Overview

Name	TNFAIP8L3 Rabbit pAb
Accession(Primary)	Q5GJ75
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q5GJ75)

Gene Name TNFAIP8L3

About Gene Symbol



Western blot analysis of TNFAIP8L3 in 293T ,
mouse brain lysates using TNFAIP8L3 antibody.

