

OAT1 Rabbit pAb

ABZ1338

Product Overview

Name	OAT1 Rabbit pAb
Accession(Primary)	O95742
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:2000; IF/ICC1:100-1:500

Target information (O95742)

Gene Name SLC22A6**About Gene Symbol**

