

SCN9A Rabbit pAb


ABZ1340

Product Overview

Name	SCN9A Rabbit pAb
Accession(Primary)	Q15858
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 226 kDa; Observed MW: 220 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (Q15858)

Gene Name	SCN9A
About Gene Symbol	



Western blot analysis of Nav1.7 in various lysates using Nav1.7 antibody.