

Smad1/5/9 Rabbit pAb


ABZ1342

Product Overview

Name	Smad1/5/9 Rabbit pAb
Accession(Primary)	Q15797
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:1000; ELISA-1:10000

Target information (Q15797)

Gene Name	SMAD1/SMAD5/SMAD9
About Gene Symbol	SMAD1 (SMAD Family Member 1) is a formin actin nucleator involved in TGF- β signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targ



Western blot analysis of Smad1/5/9 in 293T
lysates using Smad1/5/9 antibody.