

ENPP2 Rabbit pAb

ABZ1345


Product Overview

Name	ENPP2 Rabbit pAb
Accession(Primary)	Q13822
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 99 kDa; Observed MW: 99 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q13822)

Gene Name ENPP2

About Gene Symbol



Western blot analysis of ENPP2 in 293 lysates using ENPP2 antibody.

