

LAMC2 Rabbit pAb

ABZ1346

Product Overview


Name	LAMC2 Rabbit pAb
Accession(Primary)	Q13753
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 131 kDa; Observed MW: 135 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (Q13753)


Gene Name	LAMC2
About Gene Symbol	




Western blot analysis of Laminin gamma 2 in K562 lysates using Laminin gamma 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using Laminin gamma 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human lung using Laminin gamma 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human lung using Laminin gamma 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.