

IL-10 Receptor alpha Rabbit pAb

ABZ1347


Product Overview

Name	IL-10 Receptor alpha Rabbit pAb
Accession(Primary)	Q13651
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q13651)

Gene Name IL10RA

About Gene Symbol



Western blot analysis of IL-10 Receptor alpha in PC-12 lysates using IL-10 Receptor alpha antibody.

