

## PTCH1 Rabbit pAb


ABZ1348

### Product Overview


|                             |                                                                   |
|-----------------------------|-------------------------------------------------------------------|
| <b>Name</b>                 | PTCH1 Rabbit pAb                                                  |
| <b>Accession(Primary)</b>   | Q13635                                                            |
| <b>Host Species</b>         | Rabbit                                                            |
| <b>Clonality</b>            | Polyclonal                                                        |
| <b>Isotype</b>              | IgG                                                               |
| <b>Species Reactivity</b>   | Human, Mouse                                                      |
| <b>Calculated MW</b>        | Calculated MW: 161 kDa; Observed MW: 161 kDa                      |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:50-1:200; ELISA-1:10000 |

### Target information (Q13635)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | PTCH1 |
| <b>About Gene Symbol</b> |       |




Immunohistochemistry analysis of paraffin-embedded Human testis using PTCH1 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.




Immunohistochemical analysis of paraffin-embedded Human tonsils using PTCH1 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.




Immunohistochemistry analysis of paraffin-embedded Human testis using PTCH1 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Right kidney using PTCH1 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Right kidney using PTCH1 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Right kidney using PTCH1 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.