

4E BP1 Rabbit pAb


ABZ1349

Product Overview


Name	4E BP1 Rabbit pAb
Accession(Primary)	Q13541
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 13 kDa; Observed MW: 15-20 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q13541)


Gene Name	EIF4EBP1
About Gene Symbol	EIF4EBP1 (Eukaryotic Translation Initiation Factor 4E Binding Protein 1) is a protein kinase and signaling regulator involved in mTOR/PI3K/Akt signaling and translational regulation. It is an important research target in cancer, cell signaling, and therap




Western blot analysis of 4EBP1 in SH-SY5Y, C2C12 and 3T3 lysates using 4EBP1 antibody.



Immunoprecipitation analysis of 4E BP1 in HeLa lysates using 4EBP1 antibody.



Western blot analysis of 4E BP1 in Jurkat, Hela and C6 lysates using 4EBP1 antibody.



Western blot analysis of 4E BP1 in C2C12, A549 and 3T3 lysates using 4E BP1 antibody