

DUSP8 Rabbit pAb

ABZ1351

Product Overview

Name	DUSP8 Rabbit pAb
Accession(Primary)	Q13202
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB-1:1000-1:2000

Target information (Q13202)

Gene Name DUSP8**About Gene Symbol**

Western blot analysis of DUSP8 in 293T lysates using DUSP8 antibody.

