

Proteasome 26S S2/PSMD2 Rabbit pAb


ABZ1352

Product Overview

Name	Proteasome 26S S2/PSMD2 Rabbit pAb
Accession(Primary)	Q13200
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (Q13200)

Gene Name	PSMD2
About Gene Symbol	



Western blot analysis of Proteasome S2 in Jurkat
lysates using Proteasome S2 antibody