

HuB/HuD Rabbit pAb


ABZ1353

Product Overview


Name	HuB/HuD Rabbit pAb
Accession(Primary)	Q12926
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:40000

Target information (Q12926)


Gene Name	ELAVL2/ELAVL4
About Gene Symbol	




Western blot analysis of HuB/HuD in HeLa lysates using HuB/HuD antibody.



Immunofluorescence analysis of HuB/HuD in HepG2 using ELAV2/4 antibody. The picture on the right is blocked using the synthesized peptide.



Western blot analysis of HuB/HuD in Jurkat and HUVEC lysates using ELAV2/4 antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of HuB/HuD in 293 lysates
using HuB/HuD antibody