

p300 Rabbit pAb


ABZ1355

Product Overview


Name	p300 Rabbit pAb
Accession(Primary)	Q09472
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 264 kDa; Observed MW: 300 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (Q09472)

Gene Name	EP300
About Gene Symbol	



Western blot analysis of p300 in HELA, 293T, 823AV lysates using p300 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using p300 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.