

PPAR alpha Rabbit pAb


ABZ1357

Product Overview


Name	PPAR alpha Rabbit pAb
Accession(Primary)	Q07869
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200; ELISA-1:10000-1:20000

Target information (Q07869)


Gene Name	PPARA
About Gene Symbol	



Western blot analysis of PPAR alpha in 3T3 lysates using PPAR alpha antibody.



Immunohistochemistry analysis of paraffin-embedded mouse spleen tissue using PPARα antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of PPAR alpha in Hela cells using PPARα antibody (red) and Galectin3 antibody (green)