

PKC zeta Rabbit pAb


ABZ1359

Product Overview


Name	PKC zeta Rabbit pAb
Accession(Primary)	Q05513
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:5000

Target information (Q05513)


Gene Name	PRKCZ
About Gene Symbol	




Western blot analysis of PKC zeta in various lysates using PKC zeta antibody.




Western blot analysis of PKC zeta in COS7 lysates using PKC ζ antibody.



Immunofluorescence analysis of PKC zeta in NIH/3T3 cells, using PKC zeta antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded Human brain tissue using PKC zeta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of PKC zeta in COS7 lysates treated with PMA using PKC zeta antibody. The lane on the right is blocked with the synthesized peptide.