

TIE2 Rabbit pAb


ABZ1362

Product Overview

Name	TIE2 Rabbit pAb
Accession(Primary)	Q02763
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 126 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000-1:20000

Target information (Q02763)

Gene Name	TEK
About Gene Symbol	



Western blot analysis of TIE2 in CACO2 lysates
using TIE2 antibody.