

MDM2 Rabbit pAb


ABZ1363

Product Overview


Name	MDM2 Rabbit pAb
Accession(Primary)	Q00987
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:5000

Target information (Q00987)

Gene Name	MDM2
About Gene Symbol	



Western blot analysis of MDM2 in SKOV-3 lysates using MDM2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using MDM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Immunofluorescence analysis of MDM2 in A549 cells using MDM2 antibody(green).