

FBLN2 Rabbit pAb


ABZ1364

Product Overview

Name	FBLN2 Rabbit pAb
Accession(Primary)	P98095
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 127 kDa; Observed MW: 127 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:40000

Target information (P98095)

Gene Name	FBLN2
About Gene Symbol	



Western blot analysis of FBLN2 in Hela, A549 and 293T lysates using FBLN2 antibody.