

Cyclin A1 Rabbit pAb


ABZ1366

Product Overview

Name	Cyclin A1 Rabbit pAb
Accession(Primary)	P78396
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P78396)

Gene Name	CCNA1
About Gene Symbol	



Western blot analysis of Cyclin A1 in various lysates using Cyclin A1 antibody.