

## Citruline-Histone H3 ( R2+R8+R17) Rabbit pAb


ABZ1368

### Product Overview

|                             |                                               |
|-----------------------------|-----------------------------------------------|
| <b>Name</b>                 | Citruline-Histone H3 ( R2+R8+R17) Rabbit pAb  |
| <b>Accession(Primary)</b>   | P68431                                        |
| <b>Host Species</b>         | Rabbit                                        |
| <b>Clonality</b>            | Polyclonal                                    |
| <b>Isotype</b>              | IgG                                           |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                             |
| <b>Calculated MW</b>        | Calculated MW: 15 kDa; Observed MW: 15-17 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000                               |

### Target information (P68431)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | H3C1                                                                                                                                                                                                                                                            |
| <b>About Gene Symbol</b> | H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el |



Western blot analysis of Histone H3 in HeLa, 3T3 lysates using Histone H3(Citrulline R2+R8+R17) antibody.