

## PP2A alpha Rabbit pAb


ABZ1370

### Product Overview


|                             |                                                                    |
|-----------------------------|--------------------------------------------------------------------|
| <b>Name</b>                 | PP2A alpha Rabbit pAb                                              |
| <b>Accession(Primary)</b>   | P67775                                                             |
| <b>Host Species</b>         | Rabbit                                                             |
| <b>Clonality</b>            | Polyclonal                                                         |
| <b>Isotype</b>              | IgG                                                                |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                  |
| <b>Calculated MW</b>        | Calculated MW: 36 kDa; Observed MW: 36 kDa                         |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:200-1:1000; ELISA-1:5000 |

### Target information (P67775)

|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | PPP2CA |
| <b>About Gene Symbol</b> |        |



Western blot analysis of PP2A alpha in various lysates using PP2A alpha antibody.



Western blot analysis of PP2A alpha in A549 lysates using PP2A Cα antibody.

