

# Histone H4 Rabbit pAb


ABZ1371

## Product Overview


|                             |                                                                    |
|-----------------------------|--------------------------------------------------------------------|
| <b>Name</b>                 | Histone H4 Rabbit pAb                                              |
| <b>Accession(Primary)</b>   | P62805                                                             |
| <b>Host Species</b>         | Rabbit                                                             |
| <b>Clonality</b>            | Polyclonal                                                         |
| <b>Isotype</b>              | IgG                                                                |
| <b>Species Reactivity</b>   | Human, Mouse, Rat, Cat                                             |
| <b>Calculated MW</b>        | Calculated MW: 11 kDa; Observed MW: 11 kDa                         |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:200-1:1000; ELISA-1:5000 |

## Target information (P62805)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | H4C1 |
| <b>About Gene Symbol</b> |      |



Western blot analysis of Histone H4 in COS7 and HeLa lysates using Histone H4 antibody.



Western blot analysis of Histone H4 in MCF-7 lysates using Histone H4 antibody.