

RhoA Rabbit pAb


ABZ1373

Product Overview


Name	RhoA Rabbit pAb
Accession(Primary)	P61586
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (P61586)


Gene Name	RHOA
About Gene Symbol	




Western blot analysis of RhoA in HepG2 lysates using RhoA antibody.




Immunofluorescence analysis of RhoA in rat lung using Rho A antibody (red), and DAPI (blue).




Immunofluorescence analysis of RhoA in mouse lung using Rho A antibody (red), and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human stomach cancer tissue using Rho A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded rat lung tissue using RhoA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody.



Immunohistochemistry analysis of paraffin-embedded mouse lung tissue using Rho A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.