

Actin Rabbit pAb


ABZ1374

Product Overview


Name	Actin Rabbit pAb
Accession(Primary)	P60709
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P60709)


Gene Name	ACTB/POTEKP/ACTG1
About Gene Symbol	




Western blot analysis of Actin in various lysates using Actin antibody.




Western blot analysis of Actin in HeLa lysates using Actin antibody.



Western blot analysis of Actin in various lysates using Actin antibody.



Western blot analysis of Actin in various lysates using Actin antibody.



Immunohistochemistry analysis of paraffin-embedded Human skeletal muscle using Actin antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Negative control (right) obtained from antibody preabsorption by immunogen peptide.